

WOODPECKER® PROPHYLAXIS POWDER GUIDE

Old Packaging					
New Packaging					
Model		PT-S1 SUP	PT-S2 SUP+SUB	PT-S3 SUP+SUB	PT-S4 SUP+SUB
Name		FAST	SUPER	GENTLE MINI	PLUS
Composition		SODIUM BICARBONATE	ERYTHRITOL	GLYCINE	TAGATOSE
Flavor		LEMON STRAWBERRY PEACH	NATURALLY SWEET	NATURALLY SWEET	NATURALLY SWEET
Particle size		65µm & 40µm	14µm	25µm	14µm
Mohs hardness		2.5	≈2.5	2.0	2.0-2.5
Teeth	Stains and early calculus	✓	✓	⊘	✓
	Enamel	✓	✓	✓	✓
	Enamel white spots	⊘	✓	✓	✓
	Pits and fissures	⊘	✓	⊘	✓
	Dentine	⊘	✓	✓	✓
Soft tissues	Gingiva	⊘	✓	✓	✓
	Tongue and palate	⊘	✓	✓	✓
Restorations	Composite resin	⊘	✓	✓	✓
	Glass ionomer cement	⊘	✓	✓	✓
Implants	Implants and peri-implant tissues	⊘	✓	✓	✓
Orthodontics	Brackets and appliances	⊘	✓	✓	✓



Only use WOODPECKER® powders with WOODPECKER® devices. Non-original or “compatible” powders may cause damage to both the device and oral tissues.