

Hygiene & Disinfection

# YOUR DAILY ROUTINE

Enhancing your workflows  
through innovation

**riskontrol**



METALLIC TIPS

DID YOU KNOW?

It is impossible to clean correctly the inside of reusable tips due to the presence of biofilm.<sup>(1)</sup>



Metallic tip after 1 year use

- **Not effective cleaning** of reusable tips in the long-term.
- **Hygiene maintenance takes time** between two patients.



*With regard to cleaning and sterilization, disposable tips were found to be more acceptable in clinical use than metal tips.<sup>(1)</sup>*

(1) M.Inger, V.Bennani, M.Farella, F.Bennani, RD.Cannon. Efficacy of air/water syringe tip sterilization. Australian Dental Journal 2014; 59: 87-92 - R. Puttaiah et al. Rationale for using single-use disposable air/water syringe tips. Compendium. 1999;20(11):1056-1069

# riskontrol

Your must-have to perform dental care in a safe and clean environment

- **Save time & money** in your daily basis practice  
Single-use  
Avoid cleaning and sterilization steps



\*Cost assessment per dental chair for RISKONTROL® and reusable tips use for a one-year period (Considering 1 278 consultations per dental chair on a one-year basis (2016) considering that dentists are following hygiene measures with at least 1 sterilization per day  
(3) ML Silva, M. Zumpe, J. Lespinasse, M. Aulois-Griot. Cost comparison of reusable and disposable air/water syringe tips in a large French teaching hospital. Swiss Dental Journal 2017



- **Protect your patient**  
**Single-use:** reduce the risk of direct or cross-contamination  
**Minimize the spread** of viruses

- **Ensure your clinical protection**  
20 times **less risk** means **more safety and serenity** on a daily basis, for safe care and even more satisfied patients.



(2) M.Inger, V.Bennani, M.Farella, F.Bennani, RD.Cannon. Efficacy of air/water syringe tip sterilization. Australian Dental Journal 2014; 59: 87-92




your must-have



# Indispensable in your daily management

- **Match with your air/water syringe**  
RISKONTROL® is a **universal solution** thanks to its adapters.



- **Meet your expectations for your different procedures**  
**Totally dry air:** mandatory for restorative procedures

Two separate tubes guide air and water from the syringe to the end of the tip



- **Avoid risk of ingestion**  
**Safe-locking system:** patented system to protect even more your patient



- **Reach easily difficult-to-access areas and respect the patient's anatomy**  
**Easy to bend:** find the ideal curve to reach all areas  
Soft and rounded edges for the **patient's comfort:** keeping the lip away from the working area



single use

# Guaranteed daily protection

RISKONTROL® tips are suitable to all types of clinical situations

Over the past 25 years, ACTEON® has become a key player on dentistry by focusing on a mutual well-being between patient and practitioner.

RISKONTROL® tips allow to protect yourself and your patients for complete satisfaction.



“RISKONTROL® is essential for all types of treatment and is a basic element of my technical equipment. Its adaptability provides me with optimal efficacy in all clinical situations.

Dr. Gautier Dupont  
France  
Dental surgeon

# Essential in all Workflows

RISKONTROL® is the **solution** that ensures both practitioners and patients' **protection, comfort and safety** in every daily clinical procedure.



daily protection





# Discover our RISKONTROL® range

## Add color

to your patients' care experience



classic  
by **riskontrol**

- 201 664 - White
- 201 665 - Blue
- 201 666 - Yellow
- 201 667 - Green
- 201 668 - Pink

Non-sterile disposable tips (x 250)

## Adapt your tips to the taste

of your patients



art  
by **riskontrol**

- 201 700 - Mint
- 201 701 - Mandarine
- 201 702 - Liquorice
- 201 703 - Anise
- 201 704 - Black currant

Non-sterile disposable tips (x 250)

## Protect entirely your air/water syringe

thanks to the plastic sheath



total protect  
by **riskontrol**

- 201 750
- 201 751

Non-sterile disposable tips with plastic sheath (x 100)

## Choose a totally sterile workflow

workflow



sterile  
by **riskontrol**

- 201 770
- Sterile disposable tips (x 96)

Classic by RISKONTROL®, Art by RISKONTROL®, Total Protec by RISKONTROL®: Class I medical device. Sterile RISKONTROL®: Class Is medical device. CE0459 (GMED). Read the instructions on the packaging carefully before use. All the compulsory information requested for a correct use of these devices is in the Product Feature Summary available on the laboratory's website. Manufacturer: Produits Dentaires Pierre Rolland - FRANCE. Updated on: January 2023



**riskontrol**  
Essential for all clinical procedures



# For your air/water syringe


Find the adapter that matches perfectly with your equipment in our adapters guide



# universal solution

## Equip your dental office

RISKONTROL® adapters are available for more than 150 types of dental chairs.



Easy to fit

1

It's installed.

2

It's fitted.

3

It's locked.



**PRODUITS  
DENTAIRES  
PIERRE ROLLAND**

A company of ACTEON® Group

17 avenue Gustave Eiffel,  
ZI du Phare,  
33700 Mérignac, France

Tel. +33 (0) 556 340 607  
Fax. +33 (0) 556 349 292

[info@acteongroup.com](mailto:info@acteongroup.com)

[Follow-us](#)